boso medicus smart

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DE 2

**User
Instructions**

EN 36

**Mode
d'emploi**

FR 70

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IT 104

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ES 138

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Table of contents


Product features.....	38
Contents of package	39
Symbols on the blood pressure monitor.....	40
Quick guide.....	43
Introduction.....	44
 Purpose.....	45
Blood pressure values.....	46
Systolic and diastolic blood pressure value.....	46
WHO blood pressure cut-off values.....	46
Starting up.....	48
Inserting battery	48
Battery status.....	49
General instructions for self-measurement.....	50

Table of contents

Preparing to measure your blood pressure.....	52
Note on cuff sizes.....	52
Measuring your blood pressure.....	53
Attaching the cuff.....	53
Interrupting a measurement.....	57
Measurement display.....	58
Error messages.....	60
Cleaning and Disinfection.....	62
Warranty Conditions.....	63
Customer service / Disposal.....	64
Technical data.....	65
Calibration checks - Testing instructions (for trained specialist staff only).....	68

EN

Product features

EN

Contents of package

1 **boso medicus smart**
blood pressure monitor


1 LR 6 (AA) battery

1 CA 01 standard cuff

1 User
instructions

1 Case

1 Blood pressure
record card

Symbols on the blood pressure monitor

Symbol	Function/meaning
	START/STOP button
 	Serial number
	Year of manufacture
	Type BF: The device, cuff and tubing have been designed to protect the patient against electric shock.
	Reading in progress. Flashes as soon as a heartbeat was detected.
	Cardiac rhythm disorder or movement during reading.

Symbols on the blood pressure monitor

Symbol	Function/meaning
	Symbol inflation
	Symbol quick air release
	Battery fully charged.
	Battery partly discharged.
	Replace battery as soon as possible.
	Flashing! No further reading possible. Replace battery.
	Protect against liquids

EN

Symbols on the blood pressure monitor

Symbol	Function/meaning
SYS	Systolic blood pressure in mmHg
DIA	Diastolic blood pressure in mmHg
 PUL	Pulse per minute
CE 0124	Device complies with the European Medical Devices Directive.
	Device must not be discarded with household waste.
	Manufacturer
	Read instructions for use

Quick guide

- ➔ Insert battery, ensure correct positioning (P. 48)
- ➔ Press Start button  (P. 53)
- ➔ Attach cuff (P. 53)
- ➔ Inflate the cuff to approx. 180 mmHg (P. 56)

Measured systolic, diastolic and pulse values are displayed after the measurement is completed (P. 58)

This quick guide does not replace the detailed information on the use and safety of the blood pressure monitor featured in the user instructions.

Please be sure therefore to read the user instructions!

EN

Introduction

Dear Customer,
We are delighted that you have decided to purchase a bosco blood pressure monitor. The bosco brand is synonymous with optimum quality and precision and is also number 1 among professionals: 96% of all German general practitioners, physicians and internists work in practice with blood pressure instruments from bosco (API survey conducted by GfK 01/2016). The experience gained from decades of use by professionals has also been incorporated into all the patient monitors for home measurement. This instrument has also passed our strict quality control procedures and is your reliable partner for monitoring your blood pressure values.

EN

 Please read these user instructions before using the instrument for the first time, as correct blood pressure readings can only be obtained if the instrument is operated correctly. In these user instructions, „➔“ means an action by the user.

To get help concerning usage or maintenance please contact your dealer or the manufacturer (contact details can be found on the back cover of this manual).

An unexpected operating condition or an incidence which has worsened the health condition or could have worsened the health condition should be reported to the manufacturer immediately.

Introduction

If the instrument is to be sold, ensure that these user instructions are enclosed.

The bosomedicus smart blood pressure monitor is suitable for patients of all ages with an upper arm circumference of between 22 and 48 cm.

The use of accessories not referred to in these instructions for use can undermine safety.

Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 11 feet should be kept from such devices.

Purpose

Non-invasive recording of systolic and diastolic blood pressure and pulse rate in adults. Not suitable for use on newborn babies and infants. The device is suitable for self-measurement at home. The patient is regarded as the operator.

The instrument can be used with any arm circumference as shown on the corresponding cuff.

A circular logo containing the letters "EN" in a bold, sans-serif font. The logo is positioned on a dark grey rectangular background that is part of a larger design element on the right side of the page.

Blood pressure values

To determine the blood pressure, **two values** need to be measured:

- **Systolic (upper) blood pressure:**

This is produced when the heart muscle contracts and blood is forced into the blood vessels.

- **Diastolic (lower) blood pressure:**

This is present when the heart muscle expands and fills with blood again.

Blood pressure measurement values are expressed in mmHg (mm of mercury).

EN

The World Health Organization (WHO) has defined the following cut-off values for assessing blood pressure:

	systolic	diastolic
High blood pressure (hypertension)	over 140 mmHg	over 90 mmHg
Borderline normal blood pressure	130 to 139 mmHg	85 to 89 mmHg
Normal blood pressure	120 to 129 mmHg	80 to 84 mmHg
Optimum blood pressure	up to 119 mmHg	up to 79 mmHg

Blood pressure values

Whether medical treatment is necessary depends not only on the blood pressure, but also on the patient's risk profile. Please contact your GP if one of the values (systole, diastole) consistently exceeds the cut-off point for high blood pressure.

Self-measurement of blood pressure often produces somewhat lower values at home

than at the doctor's. The German Hypertension League has therefore defined lower cut-off values for home measurements:

Home measurement:	135/85 mmHg
Measurement in the doctor's surgery:	140/90 mmHg

Starting up

The device must be installed and put into operation in accordance with the instructions in this manual.

Illustration 1

- ➔ Use only a high quality leak-proof battery of the correct specifications (see "Technical data" P. 65).

⚠ If the battery is not inserted correctly, the monitor will not work and this may cause the battery to heat up and leak, damaging the monitor.

➔ Inserting the battery

The battery compartment is located on the underside of the monitor. Insert battery in the battery compartment as shown in illustration 1.

- ➔ Remove the battery if the monitor is not to be used for any length of time.

⚠ If the battery is short-circuited, it can become hot and cause burning.

Starting up

The monitor has a battery status display

-  Battery fully charged.
-  Battery partially charged.
-  Replace battery soon.
-  Flashing! No further measurements possible. Replace battery.

 Do not touch the batteries and the patient at the same time.

Care for the environment

Do not dispose of used batteries and charger units in domestic waste. You can take these to collection sites for used batteries or special waste. Contact your local authority for information.

EN

General instructions for self-measurement

EN

1. Fluctuations in blood pressure are quite normal. Even when measurements are repeated, there can be marked differences. Single or irregular measurements do not provide reliable information about the actual blood pressure. A reliable evaluation is only possible if regular measurements are taken under comparable conditions and the measurement values entered in the blood pressure record card.
2.  Self-testing does not mean self-treatment. Do not on **any** account change the dose of medicine prescribed by your doctor of your own accord.
3. Irregular heartbeats can affect the accuracy of the monitor or result in incorrect measurements (see page 61).
4. Incorrect measurements can also occur in patients with a weak pulse (possibly in patients wearing a pacemaker, for example). The blood pressure monitor itself has no effect on the pacemaker.
5. The device has not been validated for use on pregnant women.
6. Always take the measurement in a calm and relaxed state.
You are recommended to check your blood pressure twice a day – in the morning on getting up and in the evening once you have relaxed after work.
7. Always measure your blood pressure from the arm with the higher blood pressure values (where not required otherwise (see page 54)).

General instructions for self-measurement

To do this, measure the blood pressure first in both arms and thereafter always from the arm with the higher blood pressure.

Preparing to measure your blood pressure

Note on cuff sizes:

 The monitor may only be used with the following cuff types. These should be chosen according to the upper arm circumference printed on them.

Type	Arm circumference	Order number
CA01	22 – 32 cm	143-4-762*
CA02	32 – 48 cm	143-4-763

* supplied as standard.

Observe the following instructions when measuring your blood pressure:

Avoid smoking or drinking coffee for one hour before measuring your blood pressure.

Sit comfortably for the blood pressure measurement. Support your back and arms. Do not cross your legs. Place your feet flat on the floor.

Relax for 5 minutes before taking a measurement.

Do not move during the measurement.

Preparing to measure your blood pressure

- ➔ Press the Start button

All digital elements and function symbols light up momentarily to indicate that the unit is functioning correctly (see illustration 2).


Illustration 2

Attaching the cuff

The measurement should be taken from the bare upper arm.

⚠ In the case of tightly fitting outer clothing, ensure that the blood supply to the arm is not cut off when the sleeve is rolled up (if necessary remove the article of clothing before taking the measurement).

- ➔ Open out the cuff into the shape of a ring and slide it over the upper arm until the bottom edge of the cuff is about 2-3 cm above the elbow. The cuff must be positioned so that the marking is over the artery (see illustration 3).

Preparing to measure your blood pressure


Illustration 3

EN

⚠ Do not attach the cuff over wounds as this can cause further injuries.

⚠ Ensure that the cuff is not attached to an arm where arteries or veins are undergoing or have undergone medical treatment (e.g. shunt).

⚠ For women who have had a breast amputated, do not attach the cuff to the arm on the amputated side of the body.

⚠ During the measurement, medical apparatus that are being used on the same arm at the same time, may malfunction.

⚠ The cuff must not be too tight. There should be enough space to fit roughly two fingers between arm and cuff.

⚠ The metal ring should **never** lie over the artery (see illustration 3) as this may produce a false reading.

➔ Grasp the end of the cuff that feeds through the metal ring and wrap it around the outside of your arm. Attach the Velcro fastener by applying gentle pressure.

Preparing to measure your blood pressure

- ➔ Place your arm with the cuff attached on the table, holding it relaxed and slightly bent so that the cuff is level with your heart.


Illustration 4

-  Do not talk during the measurement.

-  Check that the air tube is not kinked during the measurement. This could cause congestion of the blood flow and subsequent injury.

-  The blood pressure measurement should not prevent the flow of blood for an unnecessary length of time (> 2 minutes). If the monitor fails to function correctly, remove the cuff from the arm.

-  Excessive frequent measurements can cause injury by restricting blood flow.

-  The device may not be used by unsupervised children.

-  Do not use the device near infants. This can lead to accidents or damage.

-  Do not start the device without putting on the cuff.

-  There are small parts that may cause a choking hazard if swallowed by mistake by infants.

Measuring your blood pressure

 If the ambient temperature changes significantly, leave the device about an hour to adapt to the new ambient temperature before using it.

  The performance of the device can be affected by excessive temperature, humidity or altitude.

➔ Now take the rubber ball into your right hand and inflate the cuff to approx. 180 mmHg. A triple signal tone occurs. The rising cuff pressure is shown on the digital display. In case the pressure in the cuff is not high enough, the symbol ▲ appears during the measurement. Please inflate the cuff by approx. 40 mmHg more. This procedure may be repeated until the

pressure in the cuff is higher than the systolic blood pressure (symbol ▲ disappears and the triple signal tone occurs).

At the following measurement the necessary inflation pressure is adjusted by the **memomatic** function, i.e. inflation pressure = SYSTOLIC of last measurement taken +40 mmHg). Please inflate until the triple signal tone indicated the correct pressure.

Important: The additional increase in the cuff pressure can cause faulty readings. If this occurs, repeat the measurement after a 2 minute interval and this time inflate straight to above your upper (systolic) blood pressure value.

Measuring your blood pressure

Interruption of Measurement

To interrupt a measurement you can press the button of the release valve until all air has been expelled from the cuff. (To repeat a measurement, simply reinflate the cuff.)

After inflation the air is automatically released. The decreasing cuff pressure is registered the left side of the display.

As soon as the first pulsations are registered, symbol  blinks in rhythm with the pulse.

As soon as the registered pulse is ready to be evaluated by the microprocessor the unit emits a bleep, again in pulse rhythm. At this point it is crucial that the arm with the cuff is not being moved.

The completion of the measurement is indicated by an extended bleep while the two measurement values in the display remain static.

Symbol  lights up to indicate that the quick air release can now be initiated.

➔ To do this, press the valve button until all the air has been expelled from the cuff.

EN

Measurement display

In the display appear your BP values:

SYSTOLE (upper value)

left display

DIASTOLE (lower value)

right display

Illustration 5

automatically alternating with the above display:

PULSE (per minute)

Illustration 6

If the symbol “♥” appears after the measurement, it is recommended that the measurement is repeated, holding the arm completely still. If the “♥” symbol appears again even with the arm held still, this may be an indication of an irregular heartbeat. Please discuss this with your doctor at your next visit.

Measuring your blood pressure

Blood pressure is a dynamic parameter and can be affected by the patient's position, e.g. sitting, standing, lying, moving, before or during the measurement, physical condition (stress, disease, etc.).

➔ If the measurement result is obviously wrong, repeat the measurement.

The unit switches itself off automatically after about 1 min. To repeat a measurement, press the start button  again.

After an interval of two minutes a repeat measurement can take place. Only the cuff needs to be inflated again.

➔ If no further measurement is required, remove the cuff from the arm.


Error messages

If any problems occur during the measurement, the display will show an error message instead of a blood pressure reading.

Meaning of error messages:

Err	Cause of Err message:	Troubleshooting:
	Excessive variation in cuff pressure during measurement	Hold arm still
	No valid pulse readings	Check position of cuff (see page 53) and repeat measurement.
	Systole – diastole difference too small: Systole – diastole ≤ 10 mmHg	Check position of cuff (see page 53) and repeat measurement.
Err CUF	Cause of Err CUF message:	Troubleshooting:
	Incorrectly inflated	Cuff possibly too loose, attach cuff more firmly (see page 54)

Error messages

	Cause of Err 9 message: Internal Error	Troubleshooting: Remove the batteries and briefly press START button. Put the batteries back into the monitor. If the error persists, return the instrument for repairs to the manufacturer.
	“♥” Display:	Troubleshooting: Repeat measurement recommended, holding the arm absolutely still. If the “♥” symbol appears again, even with the arm held still, this may be an indication of an irregular heartbeat. Please discuss this with your doctor at your next visit.

Cleaning and Disinfection

Cleaning Monitor:

Use a soft, dry cloth to clean your monitor.

Cleaning Cuff:

Small stains on the cuff can be removed carefully with a proprietary washing-up liquid.

Disinfection Cuff:

For disinfectant wipes (at least 5 minutes exposure time) of the device and the cuff, we recommend the disinfectant antiseptic liquid (Schülke & Mayr). To disinfect the cuff, we recommend spray disinfection. In particular, if the device is used on multiple users is to pay attention to regular cleaning and disinfections of the cuff.

Warranty Conditions

We give 3 years warranty from the date of purchase. The purchase date has to be proven by the invoice. Within the warranty period defects are eliminated free of charge. After repairs the warranty period is not extended on the whole unit, but only to the replaced components.

Excluded from the warranty are parts subject to normal wear and tear (e.g. cuff), transport damages and any damage caused by improper handling (e.g. non-compliance with the instructions for use). Damages due to disassembly by unauthorized persons are also excluded from warranty.

No claims for damages against us are substantiated by the warranty.

In the case of justified warranty claims the device has to be sent along with the original invoice to:

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Customer service

Customer service:

Warranty and repair work must be undertaken by trained and authorised personnel. Do not modify this equipment without authorization of the manufacturer.

Send the monitor, carefully packaged and with sufficient postage, to your authorised dealer or directly to:

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Disposal

Do not dispose of monitors and batteries in the domestic waste.

At the end of its lifespan, the monitor must be taken to a collection site for obsolete electronic items.

Please note that batteries and rechargeable batteries must be disposed of separately (local authority collection site).

Expected operational lifetime of the equipment: 10 years.

Expected operational lifetime of the cuff: 10.000 measurement cycles.

Technical data

Measurement principle:	oscillometric
Measurement range:	Systolic: 60 to 279 mmHg, Diastolic: 40 to 200 mmHg Pulse: 40 to 180 per minute
Cuff pressure:	0 to 299 mmHg
Display:	LCD
Operating conditions:	environmental temperature +10°C to +40°C Relative humidity 15 to 85% Air pressure 800 hPa to 1060 hPa
Transport/ Storage conditions:	environmental temperature -20°C to +60°C Relative humidity 10 to 95% Air pressure 700 hPa to 1060 hPa
Power supply:	DC 1,5 V (1 x 1.5 V Mignon IEC LR 6 alkaline manganese battery)
Typical battery life:	6.000 measurement cycles (depending on frequency of use).

Technical data

Battery check:	symbol display in LCD window
Weight:	100 g without battery
Dimensions (WxHxD):	105 mm x 35 mm x 50 mm
 Classification:	Type BF (⚡)
Protection against solid objects and water:	IP20 IP classification is the degree of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12mm diameter and greater such as fingers. This device is not protected against water.
Clinical test:	Accuracy complies with the requirements of ISO 81060 Part 2

Technical data

Maximum deviation of cuff pressure measurement: ± 3 mmHg or 2% of the reading (whichever is greater)

Maximum deviation of pulse rate display: ± 5 %

Applicable standard: IEC 80601-2-30 : "Particular requirements for basic safety and essential performance of automated non-invasive sphygmomanometers"

Calibration checks - Testing instructions

(for trained specialist staff only)

Calibration checks – every 2 years at the latest – must be performed by the following institutions or persons:

- Manufacturer
- Trade measurements authorities
- Persons who fulfil the requirements of the legislation governing the operation of medical devices.

A) Function testing

Function testing of the monitor can only be carried out on a person or with a suitable simulator.

B) Testing of Density of Pressure Circuit and Divergence from Pressure Indication

1. Remove battery.
2. Hold Start button  pressed and at the same time insert battery.
3. Release Start button . The unit now is ready for testing. The displays of SYS and DIA now show the actual pressure.
4. Prior to testing the density, open the battery compartment and take the T-connector off the unit. Wrap the cuff around a fixed cylinder and connect it with the T-connector (without deflation valve) to the pump ball and the unit.

Calibration checks - Testing instructions

(for trained specialist staff only)

Mind the setting time for cuff: at least 30 seconds.

Accomplish testing in the usual manner.

5. Prior to testing the divergence of the pressure indication, open the battery compartment and take the T-connector off the unit. Connect the unit with the normal unit and a pump ball. Accomplish testing in the usual manner.

6. Switch off unit.

7. Re-connect the cuff and the pumpball with the T-connector.

C) Placing of Safety Stamps

For complete safety cover up the drill hole for the base screw. Alternatively, connect upper housing and lower housing with a safety stamp placed directly on one of the housing edges.

 The device must not undergo service or maintenance while it is in use.

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BOSCH + SOHN GmbH u. Co. KG